


A 204
Mobile unit
for the optimum loading of Micro instruments.



- Level 1: suitable for inserts
- Level 2 right: 22 connections for hollow instruments
- Level 2 left: 11 Luer lock adapters, male, horizontal Hose
- With integrated spray arm
- Tubular filter A 800 can be used in both injector manifolds

EAN: 4002515460733 / Material number: 09862460 / Old Material Number: 69620400D

| Product category | |
|---|-----------------|
| Washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Mobile unit for ophthalmology instruments | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Option to use a central filter | • |
| Area of application | |
| Reprocessing of ophthalmology instruments | • |
| Capacity | |
| Theatre sets for micro instruments [number] | 2-4 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 365 |
| External dimensions, net width in mm | 529 |
| External dimensions, net depth in mm | 535 |
| External dimensions, gross height in mm | 530 |
| External dimensions, gross width in mm | 560 |
| External dimensions, gross depth in mm | 540 |
| Net weight in kg | 11.50 |
| Gross weight in kg | 12.50 |
| Companion pack | |
| E 476 holder [number] | 1 |
| E 477 stopper [number] | 1 |