


A 303  
Module  
for the optimum loading of up to 98 pipettes.



- Suitable for 98 measuring and full pipettes, for example
- Height of holding frame 150 mm
- 1 unit can be inserted into each of the baskets A 150 and APLW 106
- Vertical clearance without upper basket: 450 mm (with A 150 lower basket)

EAN: 4002515460825 / Material number: 09862560 / Old Material Number: 69630300D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Insert/module for pipettes	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
<b>Capacity</b>	
Pipettes [number]	98
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	185
External dimensions, net width in mm	222
External dimensions, net depth in mm	471
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	2.50
Gross weight in kg	3.50