


A 304

Module

For the optimum loading of up to 98 vials.



- Suitable for 98 centrifuge tubes, phials, test tubes, for example
- The loading heights differ depending on the positioning
- 1 unit can be inserted into each of the baskets A 100, A 150 and APLW 106

EAN: 4002515460832 / Material number: 09862570 / Old Material Number: 69630400D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
Reprocessing of test tubes	•
<b>Capacity</b>	
Phials [number]	98
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	130
External dimensions, net width in mm	222
External dimensions, net depth in mm	471
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	2.30
Gross weight in kg	3.30