

A 801
 Adapter
 for use of AUF 1 and AUF 2 in connection with A 105/1.



EAN: 4002515460931 / Material number: 09862940 / Old Material Number: 69680100D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Holder and/or adapter for transmission instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	50
External dimensions, net width in mm	50
External dimensions, net depth in mm	50
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.01
Gross weight in kg	0.01