

A 803
Holder
for the optimum loading of transmission instruments.

EAN: 4002515460993 / Material number: 09863100 / Old Material Number: 69680300D



- For use in A 105/1 upper basket along with A 800 tubular filter
- Suitable for numerous contra-angles & hand pieces as well as turbines

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Holder	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	32
External dimensions, net width in mm	32
External dimensions, net depth in mm	40
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Net weight in kg	0.05
Gross weight in kg	0.09