


KT 500

Basket carrier

For baskets with external dimensions of up to 50 x 50 cm.



- Required for PG 8172 with old Miele tables (G 720/21)
- The basket carrier is used instead of the standard basket carrier

EAN: 4002515464656 / Material number: 09876230 / Old Material Number: 67850002D

Product category	
Throughfeed tank dishwashers	•
Machine affiliation	
PG 8172	•
PG 8172 Eco	•
PTD 901	•
Product category	
Carrier/basket carrier	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, gross height in mm	345
External dimensions, gross width in mm	410
External dimensions, gross depth in mm	1270
Net weight in kg	14.00
Gross weight in kg	15.00