

A 804

Cleaning brush

For the optimum cleaning of the tubular filter A 800.



- Suitable for cleaning tubular filter A 800
- Twisted brushes with spiral trim
- End of handle with eye
- Length 500 mm, brush length 100 mm

EAN: 4002515527238 / Material number: 09960720 / Old Material Number: 68680400D

Product classification	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Filter	•
Other accessories	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	500
External dimensions, gross height in mm	45
External dimensions, gross width in mm	70
External dimensions, gross depth in mm	610
Diameter in mm	20
Net weight in kg	0.07
Gross weight in kg	0.16