

DOS G 80 ProfiLine NER

Hose dispenser pump for 5-10 l containers

With level indicator for the automatic dispensing of liquid cleaning agents.



- Suitable for all PG 8080-83 SCVi dishwashers
- 333 mm suction tube for 5 and 10 litre containers
- Dispenses up to 55 Milliliter ml per minute
- Set via the dishwasher controls
- Supplied with 4 lids for various container sizes

EAN: 4002515522936 / Material number: 09971980 / Old Material Number: 67747030NER

Product classification	
Fresh-water dishwashers	•
Machine affiliation	
PG 8130	•
PG 8130 U	•
PG 8130 i	•
PG 8130 U/i	•
PG 8131 i	•
PG 8132 SCi	•
PG 8133 SCVi	•
Product category	
External dispensing system	•
Product properties	
Material	Plastic
Colour	Grey
Production rate in ml/60 s	55
Length of suction lance in mm	333
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-230
Frequency in Hz	50
Total rated load in kW	0.01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	160
External dimensions, gross width in mm	198
External dimensions, gross depth in mm	478
Net weight in kg	2.15
Gross weight in kg	2.52