

TE  
Single use resin  
For the regeneration of partial demineralisation cartridges

EAN: 4002515551530 / Material number: 10020250 / Old Material Number: 67809510D



- 1 bag with 25 litre contents
- For empty cartridge LP 2800 & full demin. cartridge TE 2000/2800

Product classification	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 701	•
PTD 702	•

TE  
Single use resin  
For the regeneration of partial demineralisation cartridges

EAN: 4002515551530 / Material number: 10020250 / Old Material Number: 67809510D



PTD 703	•
PTD 704	•
<b>Product category</b>	
Water treatment	•
<b>Product properties</b>	
Volume in l	25
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	140
External dimensions, gross width in mm	530
External dimensions, gross depth in mm	580
Net weight in kg	18.20
Gross weight in kg	19.29