

SPP 0000

PROtect Complete

according to specific maintenance protocol incl. spare parts and repairs.

PROtect

EAN: 4002515551868 / Material number: 10077490 / Old Material Number:



- Peace-of-mind protection at fixed price
- All maintenance and repair costs covered (incl. spare parts)
- Minimised failure rates thanks to preventive maintenance
- Full cost control with a fixed contract price
- Extended appliance service life through replacement of maintenance parts

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Container and trolley washers	•
Mattress disinfectors	•
Small sterilisers	•
Large sterilisers	•
Bedpan washer-disinfectors	•
Contract contents	
Determining the technical appliance status	•
Carrying out annual maintenance work according to the checklist	•
Safety control of electrical and mechanical components	•
Carrying out safety-related measurements	•
Obtaining the necessary documentation	•
Replacing specified expendable parts	•
Spare parts for repairs	•
Travelling for selected service	•
Working time for selected service	•
Travelling for repairs	•
Working time for repairs	•
Detailed documentation	•
Electrical safety check	•
Interval	1 x annually