

PG 8596

Dispensing cabinet

For process chemicals and dispenser modules

EAN: 4002515558324 / Material number: 10087340 / Old Material Number: 69859602D



- 2 removable drawers for 3 x 5 l canisters each
- Internal dimensions H 285-690 , W 250, D 425-555 mm
- Each level 2 holding options per level for dispenser modules
- Drip tray facilitates better residual emptying of the canister

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Dispenser system	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	835
External dimensions, net width in mm	300
External dimensions, net depth in mm	600
External dimensions, gross height in mm	900
External dimensions, gross width in mm	400
External dimensions, gross depth in mm	740
Net weight in kg	26.50
Gross weight in kg	30.00