

UBS 2

Conversion kit

for use with selected series G 78 mobile units in PG 85 machines.



- For older injector units without drying unit connection

EAN: 4002515588987 / Material number: 10131100 / Old Material Number: 69215015D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PG 8583	•
PG 8593	•
Product category	
Conversion kit	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	90
External dimensions, net width in mm	129
External dimensions, net depth in mm	295
External dimensions, gross height in mm	125
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	350
Net weight in kg	0.51
Gross weight in kg	0.71