

## UBS 3

### Conversion kit

For use with selected series G 78 upper baskets in PG 85 machines.

EAN: 4002515588994 / Material number: 10131110 / Old Material Number: 69215016D



- For older load carriers without dryer connection
- For older load carriers with welded water supply

<b>Product category</b>	
Laboratory washers	•
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Conversion kit	•
<b>Product properties</b>	
Colour	Stainless Steel, Black
<b>Dimensions and weight</b>	
External dimensions, net height in mm	90
External dimensions, net width in mm	305
External dimensions, net depth in mm	176
External dimensions, gross height in mm	100
External dimensions, gross width in mm	310
External dimensions, gross depth in mm	210
Net weight in kg	0.59
Gross weight in kg	0.69