


## A 805 Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.



- For special instruments, surgical clamps (Si)

EAN: 4002515585047 / Material number: 10131950 / Old Material Number: 69680501D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Module/rack for robotic instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of robotic instruments	•
<b>Capacity</b>	
Robotics instruments [number]	1
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	92
External dimensions, net width in mm	41
External dimensions, net depth in mm	165
External dimensions, gross height in mm	25
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	485
Net weight in kg	0.37
Gross weight in kg	0.49