


A 807
Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.



- For special instruments, surgical clamps, older generation, S/Si

EAN: 4002515585078 / Material number: 10131980 / Old Material Number: 69680701D

| | |
|---|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| PG 8527 | • |
| PG 8528 | • |
| Product category | |
| Module/rack for robotic instruments | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of robotic instruments | • |
| Capacity | |
| Robotics instruments [number] | 1 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0-0 |
| Frequency in Hz | 0-0 |
| Total rated load in kW | 0.00-0.00 |
| Fuse rating in A | 0-0 |
| Dimensions and weight | |
| External dimensions, net height in mm | 77 |
| External dimensions, net width in mm | 41 |
| External dimensions, net depth in mm | 165 |
| External dimensions, gross height in mm | 25 |
| External dimensions, gross width in mm | 245 |
| External dimensions, gross depth in mm | 485 |
| Net weight in kg | 0.37 |
| Gross weight in kg | 0.49 |