

MC/2

Transfer trolley

height adjustable, for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- 2 loading racks (trough shaped) as well as handle and cover plate

EAN: 4002515588628 / Material number: 10138690 / Old Material Number: 69200109D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Storage/transport system	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in mm	1003
External dimensions, net width in mm	600
External dimensions, net depth in mm	805
External dimensions, gross height in mm	1119
External dimensions, gross width in mm	835
External dimensions, gross depth in mm	1280
Net weight in kg	34.20
Gross weight in kg	52.40