

DGGL 50-40

Perforated steam cooking containers

For DG 7000 steam ovens and DGM 6000/7000 steam ovens with microwave .

EAN: 4002515603970 / Material number: 10168190 /
Old Material Number: 23996119EU1



- For cooking vegetables, fish, meat and potatoes without liquid
- Perfect for defrosting frozen food
- Capacity: 2.2 l, usable capacity: 1.6 l
- Exterior dimensions W x D x H: 375 x 197 x 40 mm
- Not suitable for microwave operation!