

DGGL 100-40

Perforated steam cooking containers

For DG 7000 steam ovens and DGM 6000/7000 steam ovens with microwave .

EAN: 4002515604083 / Material number: 10170300 / Old Material Number: 23996124EU1



- For cooking vegetables, fish, meat and potatoes without liquid
- Perfect for defrosting frozen food
- Capacity: 4.0 l, useable capacity: 3.4 l
- Exterior dimensions W x D x H: 375 x 394 x 40 mm
- Not suitable for microwave operation!