

DOS G 80/72

Hose dispenser pump for 20-30 l containers

With level indicator for the automatic dispensing of liquid cleaning agents.



- Suitable for all PG 8172 dishwashers
- 448 mm suction tube for 20-30 litre containers
- Dispenses up to 55 Milliliter ml per minute
- Set via the dishwasher controls

EAN: 4002515648698 / Material number: 10210160 / Old Material Number: 67747022D

Product classification	
Throughfeed tank dishwashers	•
Machine affiliation	
PG 8172	•
PG 8172 Eco	•
PTD 901	•
Product category	
External dispensing system	•
Product properties	
Material	Plastic
Colour	Grey
Production rate in ml/60 s	55
Length of suction lance in mm	448
Standard electrical connection	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-230
Frequency in Hz	50
Total rated load in kW	0.01
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	160
External dimensions, gross width in mm	198
External dimensions, gross depth in mm	478
Net weight in kg	1.91
Gross weight in kg	2.13