

Set GF-1 100 mm
 Set of machine feet
 Can be adjusted to the optimum operating height.



- Suitable for all dishwashers in the PG 8172 range
- Comprising 4 machine feet
- To raise by 10 cm

EAN: 4002515641866 / Material number: 10249450 / Old Material Number: 67820030D

Product category	
Throughfeed tank dishwashers	•
Machine affiliation	
PG 8172	•
PG 8172 Eco	•
PTD 901	•
Product category	
Other accessories	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	90
External dimensions, gross width in mm	150
External dimensions, gross depth in mm	280
Net weight in kg	2.40
Gross weight in kg	2.70