

ProCare Lab 10 MA - 10 l

Liquid detergent, mildly alkaline, 10 l

for optimum reprocessing of laboratory glassware and utensils.

EAN: 4002515660300 / Material number: 10265680 / Old Material Number: 62923201EU1



- Free from surfactants and oxidants
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 11.6-11.8

Product classification	
Washer-disinfectors, laboratory	•
Applications	
Removal of blood	Good
Removal of protein	Good
Removal of residues of nutrient mediums	Good
Removal of organic residues	Good
Dispensing recommendation	
Dispensing recommendation per wash cycle	2-4 ml/l (0.2-0.4%)
pH value of solution	11.6-11.8
Main applications	
Medical laboratories	•
Biological laboratories	•
Chemical laboratories	•
Product properties	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Composition	Alkali
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for laboratory glassware	•
Limited suitability for non-ferrous metals	•
Not suitable for anodised aluminium	•
Does not contain oxidants	•
Does not contain surfactants	•
Limited suitability for non-anodised aluminium	•
Dimensions and weight	
External dimensions, net height in mm	308
External dimensions, net width in mm	191
External dimensions, net depth in mm	234
External dimensions, gross height in mm	320
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	190
Gross weight in kg	15
Fill volume in ml	10000
Compliance, certs/symbols	
GHS/CLP compliant	•