

ProCare Lab 30 P - 5 l

Neutralising agent, acidic, 5 l

For optimum neutralisation with phosphoric acid based neutralising agent.



- Free from surfactants and oxidants
- Specially developed for Miele Professional
- Suitable for all types of materials
- To be used as neutralising agent or acid pre-cleaner
- pH value 0,7

EAN: 4002515660454 / Material number: 10266490 /
Old Material Number: 62928101EU1

Product category	
Washer-disinfectors, laboratory	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing recommendation	
Dispensing recommendation per wash cycle	1-4 ml/l (0.1-0.4%)
Dispensing recommendation (per cycle)	1-4 ml/l (0.1-0.4%)
pH value at 20 °C	0,7
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Phosphate, Phosphoric acid
Dimensions and weight	
External dimensions, gross depth in mm	240.00
Product properties	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Composition	Phosphate, Phosphoric acid
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Not suitable for anodised aluminium	•
Does not contain surfactants	•
Acidic pre-cleaning	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for non-anodised aluminium	•
Not suitable for light and non-ferrous metals	•
Dimensions and weight	
External dimensions, net height in mm	217
External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	200
External dimensions, gross width in mm	140
External dimensions, gross depth in mm	240
Gross weight in kg	8
Fill volume in ml	5000

ProCare Lab 30 P - 5 l

Neutralising agent, acidic, 5 l

For optimum neutralisation with phosphoric acid based neutralising agent.

EAN: 4002515660454 / Material number: 10266490 /
Old Material Number: 62928101EU1

Compliance, test certificates, and symbols

GHS/CLP compliant

