

A 810/1

Cover net

For optimum securing of delicate items.



- Stainl. steel frames w/ plastic netting, dimensions 206 x 215 mm
- Particularly resistant and durable

EAN: 4002515671764 / Material number: 10283670 / Old Material Number: 69501800D

| Product classification | |
|---|--------------------------|
| Washer-disinfectors, medical | • |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7826 | • |
| PG 8527 | • |
| PG 8528 | • |
| G 7831 | • |
| PG 8535 | • |
| PG 8536 | • |
| PG 8562 | • |
| PG 8582 | • |
| PG 8582 CD | • |
| PG 8592 | • |
| Product category | |
| Cover | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, White |
| Dimensions and weight | |
| External dimensions, net height in mm | 0 |
| External dimensions, net width in mm | 0 |
| External dimensions, net depth in mm | 0 |
| External dimensions, gross height in mm | 20 |
| External dimensions, gross width in mm | 235 |
| External dimensions, gross depth in mm | 232 |
| Net weight in kg | 0.28 |
| Gross weight in kg | 0.35 |