

ProCare Shine 40 - 5 l
Rinsing agent, 5 l
For rinsing crockery, cutlery and glasses.



- Weakly acidic, contains citric acid
- For optimum rinsing results
- For use in commercial dishwashers
- Liquid concentrate, recommendation: 0.02-0.05 % (0.2-0.5 ml/l)
- Specially developed for Miele Professional

EAN: 4002515692691 / Material number: 10328030 / Old Material Number: 61915201EU1

Product category	
Fresh-water dishwashers	•
ProfilLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Dispensing recommendation	
Dispensing recommendation per wash cycle	0.2-0.5 ml/l (0.02-0.05%)
Large opening for suction lance	•
pH value at 20 °C	1,9 - 3,0
Recommended water range in mmol/l	0-0.53
Product properties	
Medium type	Rinse aid
Aggregate form	Liquid concentrate
Composition	Sequestering agent (no
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for aluminium	•
Suitable for silver	•
Protects decor	•
Prevents limescale and protects load items	•
Improves and shortens the drying process	•
Waste water neutral	•
Dimensions and weight	
External dimensions, net height in mm	250
External dimensions, net width in mm	150
External dimensions, net depth in mm	190
External dimensions, gross height in mm	250
External dimensions, gross width in mm	150
External dimensions, gross depth in mm	190
Gross weight in kg	5
Fill volume in ml	5000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•