

ProCare Shine 40 - 10 l

Rinsing agent, 10 l

For rinsing crockery, cutlery and glasses.

EAN: 4002515692707 / Material number: 10328040 / Old Material Number: 61915301EU1



- Weakly acidic, contains citric acid
- For optimum rinsing results
- For use in commercial dishwashers
- Liquid concentrate, recommendation: 0.02-0.05 % (0.2-0.5 ml/l)
- Specially developed for Miele Professional

Product category	
Fresh-water dishwashers	•
ProfilLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Dispensing recommendation	
Dispensing recommendation per wash cycle	0.2-0.5 ml/l (0.02-0.05%)
Large opening for suction lance	•
pH value at 20 °C	1,9 - 3,0
Recommended water range in mmol/l	0-0.53
Product properties	
Medium type	Rinse aid
Aggregate form	Liquid concentrate
Composition	Sequestering agent (no
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for aluminium	•
Suitable for silver	•
Protects decor	•
Prevents limescale and protects load items	•
Improves and shortens the drying process	•
Waste water neutral	•
Dimensions and weight	
External dimensions, net height in mm	315
External dimensions, net width in mm	190
External dimensions, net depth in mm	230
External dimensions, gross height in mm	320
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	230
Gross weight in kg	11
Fill volume in ml	10000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•