


U 409
Basket
For typical wash loads in bistros.



- For glasses, cups, plates, cake plates, saucers, etc.
- Quiver E 165 can also be used for cutlery

EAN: 4002515691168 / Material number: 10358580 / Old Material Number: 67640901D

Product category	
Tank dishwashers	•
Machine affiliation	
PG 8164	•
PTD 701	•
Product category	
Basket for mixed tableware	•
Product properties	
Material	Plastic
Colour	Blue, Grey, White
Number of rows	3
Dimensions and weight	
External dimensions, net height in mm	155
External dimensions, net width in mm	410
External dimensions, net depth in mm	410
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	1.52
Gross weight in kg	2.46