

U 403  
Basket  
for the optimum cleaning of cutlery.



- Close-meshed base
- Also suitable for sundae dishes
- Without cutlery inserts

EAN: 4002515691175 / Material number: 10358590 / Old Material Number: 67640301D

|   |                     |
|---|---------------------|
| <b>Product category</b>                 |                     |
| Tank dishwashers                        | •                   |
| <b>Machine affiliation</b>              |                     |
| PG 8164                                 | •                   |
| PTD 701                                 | •                   |
| <b>Product category</b>                 |                     |
| Basket for cutlery                      | •                   |
| Basket for saucepans/bowls, etc.        | •                   |
| <b>Product properties</b>               |                     |
| Material                                | Wire, rilsan-coated |
| Colour                                  | White               |
| <b>Dimensions and weight</b>            |                     |
| External dimensions, net height in mm   | 125                 |
| External dimensions, net width in mm    | 400                 |
| External dimensions, net depth in mm    | 400                 |
| External dimensions, gross height in mm | 170                 |
| External dimensions, gross width in mm  | 520                 |
| External dimensions, gross depth in mm  | 620                 |
| Net weight in kg                        | 2.00                |
| Gross weight in kg                      | 3.00                |