

U 407
Basket
For typical wash loads in bistros.



- For glasses, cups, plates, cake plates, saucers, etc.
- With different rail heights and distances between rails
- Quiver E 168 can also be used for cutlery

EAN: 4002515691199 / Material number: 10358620 / Old Material Number: 67640701D

Product category	
Tank dishwashers	•
Machine affiliation	
PG 8164	•
PTD 701	•
Product category	
Basket for mixed tableware	•
Product properties	
Material	Wire, rilsan-coated
Colour	White
Number of rows	5
Dimensions and weight	
External dimensions, net height in mm	92
External dimensions, net width in mm	400
External dimensions, net depth in mm	420
External dimensions, gross height in mm	170
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	620
Net weight in kg	1.94
Gross weight in kg	3.10