

E 168  
Cutlery quiver  
for the optimum loading of cutlery.



- Can be used in bistro basket U 407
- 4 compartments

EAN: 4002515691205 / Material number: 10358630 / Old Material Number: 67616801D

<b>Product category</b>	
Tank dishwashers	•
<b>Machine affiliation</b>	
PG 8164	•
PTD 701	•
<b>Product category</b>	
Insert for cutlery	•
<b>Product properties</b>	
Material	Wire, rilsan-coated
Colour	White
Number of compartments	4
<b>Dimensions and weight</b>	
External dimensions, net height in mm	102
External dimensions, net width in mm	230
External dimensions, net depth in mm	50
External dimensions, gross height in mm	150
External dimensions, gross width in mm	280
External dimensions, gross depth in mm	90
Net weight in kg	0.34
Gross weight in kg	0.70