

ProCare Shine 11 OB - 2 kg
 Powder cleaning agent, mildly alkaline, 2 kg
 For optimum cleaning of crockery, cutlery and glasses.

EAN: 4002515800508 / Material number: 10582690 / Old Material Number: 61911102EU1



- Contains active oxygen and phosphates
- Excellent cleaning performance, particularly on coffee/tea stains
- For use in commercial dishwashers
- Dispensing recommendation: 20–25 g per wash cycle or 2–5 g/l
- Specially developed for Miele Professional

Product classification	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Applications	
Removal of coffee and tea residues	Very good
Removal of starch/protein	Very good
Removal of dried on food residues	Very good
Removal of burnt-in grease	Very good
Removal of grease	Very good
Dispensing recommendation	
Dispensing recommendation per wash cycle	20–25 g / 2–5 g/l
Intensity of soiling	Strong
pH value of solution	11.0-11.5
Recommended water range in mmol/l	0-0.53
Product properties	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Powder
Composition	Alkalis,Oxygen bleach,Phosphate
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for aluminium	•
Suitable for silver	•
Protects decor	•
Protects against corrosion	•
Prevents limescale and protects load items	•
Dimensions and weight	
External dimensions, net height in mm	295
External dimensions, net width in mm	117
External dimensions, net depth in mm	117
External dimensions, gross height in mm	300
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	120
Diameter in mm	117
Weight in g	2000
Gross weight in kg	2
Compliance, certs/symbols	
GHS/CLP compliant	•