


A 207

Mobile unit

For the optimum loading of Micro instruments.



- Efficient reprocessing: fast loading and high throughput
- Top level: up to 8 trays for hollow instrument loading
- Lower level: 8 horizontal trays for further ophthal. instruments

EAN: 4002515815601 / Material number: 10632400 / Old Material Number: 69620700D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Mobile unit for ophthalmology instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Option to use a central filter	•
<b>Area of application</b>	
Reprocessing of ophthalmology instruments	•
<b>Capacity</b>	
Trays [quantity]	16
Theatre sets for micro instruments [number]	8
<b>Dimensions and weight</b>	
External dimensions, net height in mm	438
External dimensions, net width in mm	530
External dimensions, net depth in mm	542
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	4.40
Gross weight in kg	6.40
<b>Companion pack</b>	
A 311 tray holder [quantity]	2
A 800 tubular filter [quantity]	1
A 804 cleaning brush [quantity]	1