

DTD 2
Drip tray
for the secure placement of ProCare media canisters.

EAN: 4002515933695 / Material number: 10695460 / Old Material Number: 61910001D



Dimensions and weight	
External dimensions, net height in mm	62
External dimensions, net width in mm	262
External dimensions, net depth in mm	420
External dimensions, gross height in mm	74
External dimensions, gross width in mm	327
External dimensions, gross depth in mm	463
Gross weight in kg	1

- Space for 2 ProCare media canisters
- Secure place for 5 and 10 l canisters in the ProCare series
- Optimum residual emptying of the canister due to angle
- Suitable for use in an adjacent unit