

CM/1

ConductivityMeter

For measuring the conductivity level and the display for changing the cartridge.

EAN: 4002516057420 / Material number: 10889920 / Old Material Number: 69450002D



- Suitable for full demineralisation cartridges
- Optical and acoustic display for cartridge change
- 10 switch points in measuring range 0 – 199.9 µS/cm
- One-line display, connection to external LED possible
- Including wall bracket and power unit

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8060	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PTD 901	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
PG 8581	•
PG 8582	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Water treatment	•
Product properties	
Material	Plastic
Colour	Black

CM/1

ConductivityMeter

For measuring the conductivity level and the display for changing the cartridge.

EAN: 4002516057420 / Material number: 10889920 / Old Material Number: 69450002D



Dimensions and weight	
External dimensions, net height in mm	110
External dimensions, net width in mm	125
External dimensions, net depth in mm	52
External dimensions, gross height in mm	155
External dimensions, gross width in mm	310
External dimensions, gross depth in mm	350
Net weight in kg	1.56
Gross weight in kg	2.16