

A 314

Mesh tray

For the optimum reprocessing of ophthalmic lumened instruments.



- Incl. 1 injector manifold to connect ophthalmic instruments
- Mesh tray with removable lid
- Particularly suitable for A 207 mobile unit

EAN: 4002516321262 / Material number: 10979670 / Old Material Number: 69631400D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Other inserts	•
Module for ophthalmology instruments	•
<b>Area of application</b>	
Reprocessing of ophthalmology instruments	•
Reprocessing of lumened instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in mm	38
External dimensions, net width in mm	223
External dimensions, net depth in mm	151
External dimensions, gross height in mm	75
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	245
Net weight in kg	0.72
Gross weight in kg	0.85