


A 847
Injector nozzle
with support, for cleaning Erlenmeyer flasks.



- Nozzle length: 127 mm
- Nozzle Ø: 10 mm
Support Ø: 142 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516071471 / Material number: 11025600 / Old Material Number: 69684700D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	170
External dimensions, net width in mm	249
External dimensions, net depth in mm	35
External dimensions, gross height in mm	305
External dimensions, gross width in mm	300
External dimensions, gross depth in mm	200
Net weight in kg	0.12
Gross weight in kg	0.27