

A 848  
Injector nozzle  
with support, for cleaning measuring cylinders.



- Nozzle length: 350 mm
- Nozzle Ø: 10 mm  
Support Ø: 142 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516071488 / Material number: 11025610 / Old Material Number: 69684800D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	350
External dimensions, net width in mm	20
External dimensions, net depth in mm	20
External dimensions, gross height in mm	350
External dimensions, gross width in mm	180
External dimensions, gross depth in mm	50
Net weight in kg	0.30
Gross weight in kg	0.53