

A 306/1

Module

for the optimum loading of large-volume laboratory glassware.



- Suitable for 4 1–2 litre measuring cylinders, high shape
- Suitable for 2 1–2 litre measuring cyl., low shape & 2 high shape
- Plastic-coated contact surfaces
- Measuring cylinders even removable when module is in machine
- Maximum length of measuring cylinders: 500 mm
- Up to 2 modules can be inserted in the baskets A 150 and APLW 106

EAN: 4002516076254 / Material number: 11035210 / Old Material Number: 69630601D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Insert/module for measuring cylinder	•
Module for large-volume laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black, White
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
<b>Capacity</b>	
Measuring cylinders [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in mm	235
External dimensions, net width in mm	418
External dimensions, net depth in mm	471
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	2.10
Gross weight in kg	3.10