

## ProCare Shine 12 GC - 80 Tabs

80 detergent tabs

For optimum cleaning of crockery, cutlery and glasses.

EAN: 4002516098409 / Material number: 11037820 / Old Material Number: 61914201EU1



- Particularly simple handling thanks to water-soluble wrapping
- Contains active oxygen
- Free from phosphates and chlorine
- For light to medium soiling
- For use in commercial dishwashers

Product category	
Fresh-water dishwashers	•
ProfilLine dishwashers	•
Applications	
Removal of coffee and tea residues	Good
Removal of starch/protein	Good
Removal of dried on food residues	Good
Removal of burnt-in grease	Good
Removal of grease	Very good
Dispensing recommendation	
Dispensing recommendation per wash cycle	1 Tab
Intensity of soiling	Light to moderate
pH value of solution	10.2
Recommended water range in mmol/l	0-0.53
Product properties	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Tab
Water-soluble film	•
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for aluminium	•
Suitable for silver	•
Protects decor	•
Protects against corrosion	•
Prevents limescale and protects load items	•
Does not contain phosphates	•
Dimensions and weight	
External dimensions, net height in mm	230
External dimensions, net width in mm	190
External dimensions, net depth in mm	93
External dimensions, gross height in mm	230
External dimensions, gross width in mm	90
External dimensions, gross depth in mm	190
Number of detergent tabs	80
Weight in g	1280
Gross weight in kg	1
Compliance, test certificates, and symbols	
GHS/CLP compliant	•
EU eco label	•
Verification numbers	
EU Eco label: AT/038/006	•