


APFD 201

Use

for the optimum loading of breathing masks.



- Suitable for holding 6 breathing masks securely
- Accessories: 1 small-part basket with 6 coloured compartments
- With separate lid for particularly delicate parts (6 compartments)

EAN: 4002516080480 / Material number: 11045760 / Old Material Number: 67020101D

|   |      |
|---|------|
| <b>Product category</b>                 |      |
| Fresh-water dishwashers                 | •    |
| <b>Machine affiliation</b>              |      |
| PG 8063                                 | •    |
| <b>Product category</b>                 |      |
| Insert for breathing masks              | •    |
| <b>Area of application</b>              |      |
| Reprocessing of breathing masks         | •    |
| <b>Capacity</b>                         |      |
| Breathing masks [number]                | 6    |
| <b>Dimensions and weight</b>            |      |
| External dimensions, net height in mm   | 219  |
| External dimensions, net width in mm    | 470  |
| External dimensions, net depth in mm    | 480  |
| External dimensions, gross height in mm | 255  |
| External dimensions, gross width in mm  | 525  |
| External dimensions, gross depth in mm  | 570  |
| Net weight in kg                        | 5.70 |
| Gross weight in kg                      | 7.11 |