

A 312  
Module  
For reprocessing viscometers.



- Suitable for up to 8.0 viscometers
- Equipped with 8 x A 867 special nozzles
- Nozzle spacing 118 mm across the width, 118 mm across the depth

EAN: 4002516087397 / Material number: 11054470 / Old Material Number: 69631200D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Module for viscometers	•
<b>Product properties</b>	
Material	Stainless steel, Plastic, Silicone
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Viscometer [quantity]	8.0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	66
External dimensions, net width in mm	133
External dimensions, net depth in mm	479
External dimensions, gross height in mm	250
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	505
Net weight in kg	1.58
Gross weight in kg	2.33