

A 313

Module

For reprocessing large-volume laboratory glassware.



- Laboratory flasks, up to 10 l
- Round flasks, up to 6 l
- Erlenmeyer flasks, up to 5 l
- Measuring flasks, up to 5 l
- Can be equipped with nozzles & holders based on individual needs
- Equipped with 5 injector positions for 2 or 3 load items

EAN: 4002516087854 / Material number: 11055970 / Old Material Number: 69631300D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
Product category	
Insert/module for narrow-necked glassware	•
Module for large-volume laboratory glassware	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	221
External dimensions, net width in mm	234
External dimensions, net depth in mm	495
External dimensions, gross height in mm	260
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	505
Net weight in kg	2.82
Gross weight in kg	3.29