

## A 315 Module

For increased capacity and reprocessing of lumened instruments.



- Increases the injector capacity from 22 to 44
- Backwards-compatible with all A 151 products on the market
- Compatible with all standard dental & ENT accessories

EAN: 4002516087861 / Material number: 11055980 / Old Material Number: 69631500D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert for ENT instruments	•
Insert/module for dental instruments	•
Other inserts/other modules	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Option to use a central filter	•
<b>Area of application</b>	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
<b>Capacity</b>	
ENT instruments [number]	22
Transmission instruments [number]	22
Siphons [number]	22
Lumened instruments [number]	22
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	145
External dimensions, net width in mm	492
External dimensions, net depth in mm	542
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	1.35
Gross weight in kg	2.35
<b>Companion pack</b>	
A 838 set of spanners [number]	1