

A 300/2

Injector module

for optimum loading of laboratory glassware – can be equipped with 8 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 8 injector nozzles
- Part of the EasyLoad system
- Nozzle spacing 118 mm across the width, 118 mm across the depth
- Particularly suitable for A 840 and A 841 (nozzle diameter 6 mm)

EAN: 4002516087632 / Material number: 11056770 / Old Material Number: 69630002D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	8
<b>Standard electrical connection</b>	
Number of phases	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	72
External dimensions, net width in mm	133
External dimensions, net depth in mm	479
External dimensions, gross height in mm	200
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	0.56
Gross weight in kg	1.64