

A 860  
Holder rack  
for A 300/2 and A 300/3 injector modules.



- Arrangement: 2 x 4 fields
- For use on A 300/2 and A 300/3 injector modules
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516087663 / Material number: 11056800 / Old Material Number: 69686000D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	8
Compartment width in mm	114
Compartment length in mm	114
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	30
External dimensions, net width in mm	250
External dimensions, net depth in mm	350
External dimensions, gross height in mm	35
External dimensions, gross width in mm	300
External dimensions, gross depth in mm	400
Net weight in kg	0.10
Gross weight in kg	0.18