


APLW 002

Mobile unit

for optimum loading of large-volume laboratory glassware, 5 and 10 l.

EAN: 4002516115304 / Material number: 11104490 / Old Material Number: 69300200D



- Equipped with 5 nozzles
- For holding up to 4 flasks, max. Ø 240 mm
- For holding up to 5 flasks, max. Ø 190 mm

| | |
|--|--------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Module for large-volume laboratory glassware | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Flasks [number] | 5 |
| Dimensions and weight | |
| External dimensions, net height in mm | 413 |
| External dimensions, net width in mm | 535 |
| External dimensions, net depth in mm | 582 |
| External dimensions, gross height in mm | 580 |
| External dimensions, gross width in mm | 565 |
| External dimensions, gross depth in mm | 615 |
| Net weight in kg | 5.94 |
| Gross weight in kg | 8.49 |