

APLW 005  
 Mobile unit  
 for optimum loading of up to 56 pipettes.



- For holding up to 56 pipettes arranged diagonally

EAN: 4002516115434 / Material number: 11104520 / Old Material Number: 69300500D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Mobile unit for pipettes	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
<b>Capacity</b>	
Pipettes [number]	56
<b>Dimensions and weight</b>	
External dimensions, net height in mm	308
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	570
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	6.40
Gross weight in kg	9.20