

APLW 007  
Mobile unit  
for optimum loading of 2 pipette cartridges.



- For lower level
- Maximum pipette length: 520 mm

EAN: 4002516115458 / Material number: 11104540 / Old Material Number: 69300700D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6111	•
<b>Product category</b>	
Mobile unit for pipettes	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	580
External dimensions, net width in mm	535
External dimensions, net depth in mm	420
External dimensions, gross height in mm	580
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	610
Net weight in kg	7.18
Gross weight in kg	9.98