

APLW 011

Upper basket

with 110 unoccupied injector slots for nozzles with Ø 2.5 mm.



**Image
in progress**

- Only suitable for use on upper levels

EAN: 4002516115496 / Material number: 11104580 / Old Material Number: 69301100D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Phials [number]	110
Narrow-necked glassware [number]	110
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	532
External dimensions, net depth in mm	560
External dimensions, gross height in mm	350
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	610
Net weight in kg	5.20
Gross weight in kg	7.80