

APLW 013/2

Upper basket

equipped with 32 APLW 051 injector nozzles and 10 APLW 055 injector nozzles



**Image  
in progress**

- Only suitable for use on upper levels

EAN: 4002516115533 / Material number: 11104620 / Old Material Number: 69301302D

<b>Product category</b>	
Large-chamber laboratory washers	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Capacity</b>	
Narrow-necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
<b>Dimensions and weight</b>	
External dimensions, net height in mm	250
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	355
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	615
Net weight in kg	4.55
Gross weight in kg	9.01
<b>Companion pack</b>	
APLW 051, injector nozzle Ø 4 mm, length 110 mm [number]	32
APLW 055, injector nozzle Ø 6 mm, length 175 mm [number]	10