

APLW 014
Upper basket
with 20 unoccupied injector slots.



- Only suitable for use on upper levels

EAN: 4002516115540 / Material number: 11104630 / Old Material Number: 69301400D

| | |
|--|------------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Upper basket | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Blue, Stainless Steel, White |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow-necked glassware [number] | 20 |
| Flasks [number] | 20 |
| Erlenmeyer flasks [number] | 20 |
| Round flasks [number] | 20 |
| Dimensions and weight | |
| External dimensions, net height in mm | 21 |
| External dimensions, net width in mm | 20 |
| External dimensions, net depth in mm | 16 |
| External dimensions, gross height in mm | 340 |
| External dimensions, gross width in mm | 555 |
| External dimensions, gross depth in mm | 615 |
| Net weight in kg | 4.22 |
| Gross weight in kg | 6.65 |